

Urinometer model CLINITEC 500 SIEMENS-GJERMANI



Pershkrim:

Shjegim i sistemit Testet e matura	Aparat per matjen e urines Gjysem automatizuar Leukocite, Nitrite, Proteine, Gjak, Glukoze, Ketone, Bilirubin, Urobilinogen, pH, Pesha specifike, kreatinine, and (Protein-to- Creatinine Ratio)
Matje automatike Forma e testit Matja e testit	Ngjyren e urines Strip reagent i thate Matja e ndryshimit te ngjyres ne menyre fotometrike Lexim i dyfishte ne baze te gjatesise se vales Rregullimi automatic i ngjyres se urines

Te pergjithshme

Permasat	12.75 (h) x 15.75 (w) x 13.75 (d) in. 32 (h) x 39 (w) x 35 (d) cm
Pesha	7.2 kg (16 lbs)
Fuqia e kerkuar	100 V to 240 V; 50/60 Hz
Rryma	<0.5 mA in normal condition <3.5 mA in single fault condition
Fuqia maksimale ne hyrje	72 VA
Temperatura e ambientit te punes	18°C to 30°C (64°F to 86°F)
Lageshtia e ambientit te punes	20% to 80% relative humidity

Kalibrimi

Kalibrimi	Kalibrim automatik ne cdo test
-----------	--------------------------------

Kompujetri i brendeshem

Memorja	Ruajtja e 500 pacienteve plus 200 kontrolleve
Nderfaqja	ASTM (CCS)
Ekрани	Color LCD touchscreen, 5.7 in (14.5 cm), 320 horizontal by 240 vertical pixels (QVGA)